ARTERIOCYTE

Helping Patients Heal Faster™

MAGELLAN®

Aspiration kit

Aspirationssæt

Aspirationskit

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Australian Sponsor

Emergo Australia
Level 20 Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000, Australia

© Arteriocyte Medical Systems, Inc.
All Rights Reserved

AMS95300LBL Rev. 3

EXPLANATION OF SYMBOLS / SYMBOLFÖRKLARING / ERLÄUTERUNG DER
SYMBOLS / ΕΠΕΞΗΓΗΣΗ ΣΥΜΒΟΛΩΝ / EXPLICACIÓN DE SÍMBOLOS /
SYMBOLIEN SELITYS / EXPLICATION DES SYMBOLES / SPIEGAZIONE DEI
SIMBOLI / UITLEG VAN DE SYMBOLEN / EXPLICAÇÃO DOS SÍMBOLOS /
FÖRKLARING AV SYMBOLER / SEMBOLLERIN AÇIKLAMASI

Caution, Consult Accompanying Documents / Forsigtig, se de
medfølgende documenter / Achtung, Begleiddokumente lesen /
Προσοχή, Συμβουλευτείτε τα συνοδευτικά έγγραφα /
Precaución. Consultar la documentación adjunta / Huomio,
tutustu liiteasiakirjoihin / Attention, consulter la documentation
jointe / Attenzione consultare i documenti di accompagnamento
/ Let op, raadpleeg de meegeleverde documenten / Cuidado!
Consultar os documentos anexos / Varning, se medföljande
dokument / Dikkat, Birlikte Verilen Belgelere Başvurun

Lot Number / Partinummer / Chargennummer / Αριθμός
παρτίδας / Número de lote / Eränumero / Numéro de lot /
Numero di lotto / Partijnummer / Número de lote /
Partinummer / Parti Numarası

Catalog Number / Katalognummer / Katalognummer /
Αριθμός καταλόγου / Número de catálogo / Luettelonumero /
Référence / Numero di catalogo / Catalogusnr / Número de
catálogo / Katalognummer / Katalog Numarası

Sterilized Using Ethylene Oxide / Steriliseret ved anvendelse
af ethylenoxid / Mit Ethylenoxid sterilisiert / Αποστειρωμένο με
χρήση αιθυλενοξειδίου / Esterilizado con óxido de etileno /
Steriloitu käyttäen etyleenioksidia / Stérilisé à l'oxyde
d'éthylène / Sterilizzato mediante ossido di etilene /
Gesteriliseerd met behulp van ethyleenoxyde / Esterilizado
com óxido de etileno / Steriliserad med etylenoxid / Etilen
Oksit Kullanılarak Sterilize Edilmiştir

Latex Free / Latexfri / Latexfrei / Δεν περιέχει λάτεξ / Sin
latex / Lateksiton / Sans latex / Privo di lattice / Bevat geen
latex / Sem látex / Latexfri / Lateks İçermez

Nonpyrogenic / Ikke pyrogen / Nicht pyrogen / Μη
πυρετογόνο / Apirógeno /
Ei-pyrogeneinen / Apyrogène / Non pirogeno / Pyrogeenvrij /
Não pirogénico / Pyrogenfri / Pirojenik Değildir

Do Not Reuse / Må ikke genbruges / Nicht wiederverwenden /
Μην επαναχρησιμοποιείτε / No reutilizar / Ei saa käyttää
uudelleen / Ne pas réutiliser / Non riutilizzare / Niet
hergebruiken / Não reutilizar / Får ej återanvändas / Tekrar
Kullanmayın

Use By / Anvendes inden / Haltbar bis / Ημερομηνία λήξης /
Fecha de caducidad / Käytettävä viimeistään / Date limite
d'utilisation / Utilizzare entro il / Uiterste gebruiksdatum / Data
limite de utilização / Används före / Son Kullanma Tarihi

Open Here / Åbn her / Hier öffnen / Ανοίξτε εδώ / Abrir aquí /
Avaa tästä / Ouvrir ici / Aprire qui / Hier openen / Abrir aquí /
Öppna här / Buradan Açın

Prescription only / Receptpligtig / Verschreibungspflichtig /
Μόνο με ιατρική συνταγή / Solo con receta / Vain
hoitomääräyksellä / Uniquement sur ordonnance / Solo su
prescrizione / Uitsluitend op recept / Exclusivamente sob
receita médica / Endast efter receptförskrivning / Reçete ile
satılır

CE marking / CE-mærkning / CE-Kennzeichnung / Σήμανση
CE / Marcado CE / CE-merkintä / Marquage CE /
Contrassegno CE / CE-keurmerk / Marcação CE / CE-
märkning / CE işareti

Authorized European representative / Autoriseret europæisk
repræsentant / Autorisierter Europa-Vertreter /
Εξουσιοδοτημένος αντιπρόσωπος για την Ευρώπη /
Representante europeo autorizado / Valtuutettu edustaja
Euroopassa / Représentant européen agréé / Rappresentante
europeo autorizzato / Erkende Europese vertegenwoordiger /
Representante autorizado na Europa / Auktoriserad europeisk
representant / Yetkili Avrupa temsilcisi

Magnetic Resonance (MR) unsafe / Usikker ved magnetisk
resonans (MR) / Nicht sicher für Kernspintomographie (MR) /
Μη ασφαλές για μαγνητική τομογραφία (MR) / No seguro para
resonancia magnética (RM) / Ei turvallinen
magneettikuvauksessa / Incompatible avec la résonance
magnétique (RM) / Non sicuro per la risonanza magnetica
(MR) / Niet veilig bij magnetische resonantie (MRI-scan) /
Inseguro para Ressonância Magnética (RM) / Ej säker vid
magnetresonans (MR) / Manyetik Rezonans (MR) için güvenli
değildir

Consult instructions for use / Se brugsanvisningen /
Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες
χρήσης / Consultar las instrucciones de uso / Perehdy
käyttöohjeisiin / Consulter le mode d'emploi / Consultare le
istruzioni per l'uso / Raadpleeg de gebruiksaanwijzing /
Consultar instruções de utilização / Läs bruksanvisningen /
Kullanım kılavuzuna başvurun

Manufacturer / Producent / Hersteller / Κατασκευαστής /
Fabricante / Valmistaja / Fabricant / Produttore / Fabrikant /
Fabricante / Tillverkare / Üretici

This product contains the Phthalate, DEHP (Di(2-ethyhexyl) phthalate) in component AMS40005 / Dette produkt indeholder phthalat, DEHP (Di(2-ethyhexyl) phthalat) i komponent AMS40005 / In Komponente AMS40005 dieses Produkts ist das Phthalat DEHP (Di(2-ethyhexyl) phthalat) enthalten / Αυτό το προϊόν περιέχει την φθαλική ουσία DEHP [φθαλικό δι(2-αιθυλεξυλ) εστέρα] στο συστατικό μέρος AMS40005 / Este producto contiene el ftalato, DEHP o di(2-etilhexil) ftalato, en el componente AMS40005 / Tämä tuote sisältää ftalaattia, DEHP (di(2-etyyliheksyyli)ftalaattia) osassa AMS40005 / Ce produit contient un phtalate, le DEHP (di(2-éthylhexyl) phtalate) dans le composant AMS40005 / Questo prodotto contiene ftalato, DEHP (ftalato di bis[2-etilesile]), nel componente AMS40005 / Dit product bevat het ftalaat DEHP (bis-(2-ethyhexyl)-ftalaat) in onderdeel AMS40005 / Este produto contém ftalato, DEHP (di(2-etil-hexil)ftalato) no componente AMS40005 / Denna produkt innehåller ftalat, DEHP (di(2-etytexyl)-ftalat) i komponent AMS40005 / Bu ürün, AMS40005 bileşeninde Ftalat, DEHP (Di(2-etileksil) ftalat) içermektedir

TABLE OF CONTENTS

MAGELLAN® Aspiration Kit.....	6
INTRODUCTION	6
Indications for Use	6
Contraindications for Use	6
Cautions.....	6
Precautions	6
Warnings.....	7
How Supplied.....	7
200 Micron Blood Component Filter	7
INSTRUCTIONS FOR USE	7
Intended Use.....	7
Directions	7
Lightning Needle™ Bone Marrow Aspiration Set	8
INSTRUCTIONS FOR USE	8
Intended Use.....	8
Directions	9

INSTRUCTIONS FOR USE

MAGELLAN® Aspiration Kit

INTRODUCTION

Indications for Use

The MAGELLAN® Aspiration Kit is intended for use only with the MAGELLAN® Autologous Platelet Separator Instrument.

Note: Refer to the MAGELLAN® Autologous Platelet Separator Operator's Manual for Usage Instructions, Disclaimers, Warnings, Precautions, and Contraindications.

Contraindications for Use

- Fracture
- Excessive Tissue (Severe obesity and/or absence of adequate anatomical landmarks)
- Osteoporosis
- Infection at the area of insertion
- Previous, significant orthopedic procedure at the site, prosthetic limb or joint
- Intraosseous at site in the past 48 hours

Cautions

- Federal Law (USA) restricts this device to be sold by or on the order of a physician.
- Refer to the System Operator's Manual supplied with the instrument, as well as Instructions for Use for components from other manufacturers inside this convenience kit before performing the procedure.
- Use only MAGELLAN® disposable accessories.
- Inspect the kit prior to use. Do not use the kit if any component of the tray is damaged or opened.
- Treat all blood and fluids using universal precautions for bloodborne pathogens.
- All components of this kit are sterile fluid pathway single patient use. Do not re-sterilize.

Precautions

- Only Arteriocyte Medical Systems, Inc. sterilized disposable kits are approved for use with the MAGELLAN® Autologous Platelet Separator Instrument. It is important that aseptic technique be used to minimize the possibility of contamination of the disposable components and/or patient.
- Store all disposable components in a dry place away from extremes of environmental conditions.
- Materials used in the MAGELLAN® Disposable Kit may be sensitive to chemicals (such as solvents and certain detergents).
- Under certain adverse conditions, exposure to these chemicals (including vapors) may cause the plastics to fail or malfunction. Visually inspect the contents of the disposable kit. Should any evidence of damage to the components be found during inspection or setup, do not use the disposable components. Do not use silicone oils or greases near disposable components.

Warnings

- Check skin thickness on all patients prior to manual Cannula Assembly insertion.
- Special caution must be exercised for patients with BMI greater than 30.
- Reuse of this single-patient-use device creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

How Supplied

Each MAGELLAN® Aspiration Kit contains the following components:

- One (1) Set of Lightning Bone Marrow Aspiration Cannula and Plunger
- One (1) Lightning needle Adapter
- One (1) 60-mL Syringe
- Two (2) 30-mL Syringes
- One (1) 200 Micron Blood Component Filter
- One (1) 3-Way Stopcock
- Four (4) Sterile Syringe Caps
- One (1) 18-Gauge Needle

200 Micron Blood Component Filter

INSTRUCTIONS FOR USE

Intended Use

To filter clots and other particles from blood and blood components for delivery.

Directions

IMPORTANT: Use aseptic technique. Single patient use.

1. The filter should be transferred to the sterile field using proper technique.
2. Remove the filter and close the clamp on the filter.
3. Remove the vent cap from the distal (bottom) end of the filter (clamp end).
4. Luer connect the male port of a three-way stopcock to the open, distal end of the filter.
5. Luer connect a 60mL syringe to the female port of the three-way stopcock that is directly across (180°) from the port attached to the filter.
6. Luer connect a syringe (10mL minimum) to the side female port of the three-way stopcock. Label this syringe with "Heparin". (Heparin not included in the kit)
7. Draw 10mL of Heparin (concentration of 1000 units/mL) into an additional Heparin-labeled syringe.

IMPORTANT: Ensure that the patient record has been checked for history of Heparin allergy. For patients with a Heparin allergy, ACD-A should be used as the anticoagulant in place of Heparin.

8. Remove the proximal vent cap from the filter assembly and luer connect the Heparin-filled syringe to the open, proximal end of the filter assembly.

INSTRUCTIONS FOR USE

9. Holding the filter upright (three-way stopcock on the bottom), push the Heparin into the filter assembly and rock it back and forth to coat well, leaving Heparin in place until just before filtration of blood component.

NOTE: If filtration is being done off of the sterile field, the filter assembly should be transferred to the processing technician at this time.

10. Rotate the three-way stopcock so that the off position is pointing towards the 60mL syringe and open the distal tubing clamp.
11. Holding the filter upright (three-way stopcock on the bottom), remove the Heparin from the filter by pulling it into the Heparin-labeled syringe attached to the side port of the three-way stopcock.
12. Rotate the three-way stopcock so that the off position is pointing towards the side port.
13. Remove the empty Heparin-labeled syringe from the proximal (top) end of the filter and connect the blood component syringe to the open proximal end of filter.
14. Holding the filter upright (three-way stopcock on the bottom), simultaneously push the blood component syringe and pull the distal 60mL syringe gently to draw the blood component through the filter assembly and into the distal 60mL syringe.

NOTE: Ensure a steady flow of blood component from filter into distal tubing to avoid air from entering the tubing.

15. Once filtration is complete, rotate the three-way stopcock so that the off position is 45° between two of the ports (closing all three ports) and close the clamp on the distal (bottom) end of the filter.
16. While pulling back gently on the plunger of the 60mL filtered blood component syringe, disconnect the syringe from the filter.
17. Cap the 60mL syringe containing the filtered blood component with a sterile syringe cap.


Manufactured by: Arterocyte Medical Systems, Inc., 45 South Street, Hopkinton, MA 01748 USA

Lightning Needle™ Bone Marrow Aspiration Set

INSTRUCTIONS FOR USE

Intended Use

The Arterocyte Medical Systems Lightning™ Bone Marrow Aspiration Cannula System is intended for use in the intraosseous aspiration of viscous bone marrow material from the human body.

Directions

1. The Cannula Assembly, Plunger and Adapter should be transferred to the sterile field using proper technique.
2. Palpate and locate the entry point.
3. Using aseptic technique, prep the skin. If desired, apply local anesthetic to the prepped area.
4. Remove the Lightning Needle Cannula Assembly (Cannula and Trocar), Plunger, and Driver Adaptor from sterile package(s).
IMPORTANT: Needle is sharp, do not remove protective cover.
5. Collect 10mL of Heparin (concentration of 1000 units/mL) from a Heparin-labeled cup on the sterile field into a Heparin-labeled luer syringe.
IMPORTANT: Ensure that the patient record has been checked for history of Heparin allergy. For patients with a Heparin allergy, ACD-A should be used as the anticoagulant in place of Heparin.
6. Using caution to avoid the sharp end of the Trocar, aseptically remove (unthread) the Trocar from Cannula, leaving the Trocar in the sterile field for reassembly (see step 9).
7. Connect the Heparin-containing syringe to the threaded end of the Cannula. Flush Cannula with 10mL of Heparin while dispensing the Heparin back into the Heparin-labeled cup.
8. Once flush is complete, disconnect syringe from Cannula.
9. Using caution to avoid the sharp end of the Trocar, thread the Trocar back into Cannula ensuring the groves (threading) on Trocar hub align with groves on Cannula hub.
10. Securely attach the Driver Adapter to the Surgical Drill.
11. Secure Cannula Assembly to the Driver Adapter on the Surgical Drill.
IMPORTANT: Be careful to maintain sterile technique.
IMPORTANT: Control the patient's movement prior to and during Cannula Assembly insertion.
12. Carefully remove the protective cover from the Cannula Assembly.
13. Position Cannula Assembly over the entry point at the appropriate angle to the bone.
14. Slightly power or manually press Cannula Assembly tip through the skin until the tip touches bone.
15. Activate the Surgical Drill and insert the Cannula Assembly into the bone with gentle, steady pressure until the desired depth is attained.
16. While applying gentle downward pressure to the Cannula Assembly, remove the Surgical Drill and Adapter from the Trocar hub and set aside.
17. While applying gentle downward pressure to the Cannula, carefully remove (unthread) the Trocar from the Cannula by turning the Trocar hub one and a half turns counterclockwise and pulling back gently. Confirm Cannula stability.
18. Insert the Plunger into the Cannula and seat fully.
19. While applying gentle downward pressure to the Cannula, gently remove the Plunger from the Cannula.
20. Thread a Heparin-coated syringe onto the Cannula hub and while applying gentle downward pressure to the cannula, begin marrow aspiration by gently pulling back on the syringe plunger.
IMPORTANT: Do not pull back on the syringe body as this may destabilize the Cannula.

INSTRUCTIONS FOR USE

21. Repeat step 20 as needed to collect the amount of bone marrow desired.
22. Once marrow aspiration collection is complete, remove the collection syringe (unthread) and cap the syringe with a sterile syringe cap.
23. Attach another syringe to the Cannula and remove the Cannula from the patient by gently applying outward traction.
IMPORTANT: Do not rock or bend the Cannula during removal.
24. Discard Cannula, Trocar, Plunger and Adaptor into the proper sharps container once utilization is complete.
25. Dress insertion site as appropriate.

Manufactured by: Arteriocyte Medical Systems, Inc., 45 South Street, Hopkinton, MA 01748 USA

DISPOSABLES LIMITED WARRANTY – ARTERIOCYTE MEDICAL SYSTEMS, INC. (OUTSIDE U.S.)¹

THE FOLLOWING LIMITED WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES.

- A. This LIMITED WARRANTY provides assurance for the customer who receives an Arteriocyte Medical Systems® ("AMS") **Magellan® Aspiration Kit** ("Product"), that should the Product fail to function to specification, AMS will issue a credit equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any AMS replacement Product used for that customer. THE WARNINGS CONTAINED IN THE PRODUCT LABELLING ARE CONSIDERED AN INTEGRAL PART OF THIS LIMITED WARRANTY. CONTACT YOUR LOCAL AMS REPRESENTATIVE TO OBTAIN INFORMATION ON HOW TO PROCESS A CLAIM UNDER THIS WARRANTY.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
 - (1) The Product must be used prior to its "Use By" date.
 - (2) The Product must be returned to AMS within sixty (60) days after use and shall be the property of AMS.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must not have been used more than one time for any customer.
 - (5) The Product must be used in conformity with the user instructions in the "User Manual" delivered to the customer together with the Product, of which this LIMITED WARRANTY is an integral part.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
 - (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper use or material alteration of the replaced Product.
 - (2) AMS is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort, patent infringement or otherwise.
- D. This LIMITED WARRANTY does not cover those parts that, by their very nature, are likely to deteriorate or which AMS considers should be periodically replaced consistently with normal routine or preventative maintenance requirements.
- E. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part to term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.
- F. No representative, agent, dealer, retailer, or intermediary of AMS shall have authorization to amend the contents of this LIMITED WARRANTY.
- G. The validity, interpretation and performance of the agreement for which this LIMITED WARRANTY is issued, as well as any disputes relating or connected thereto is controlled by and construed under the laws of the State of Delaware, USA.

¹This Limited Warranty is provided by the Manufacturer, Arteriocyte Medical Systems, Inc., 45 South Street, Hopkinton, MA, 01748, USA. It applies only outside of the United States. Areas outside the United States should contact their local Arteriocyte Medical Systems representative for exact terms of the Limited Warranty

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Üretici: Arterioocyte Medical Systems, Inc., 45 South Street, Hopkinton, MA 01748 ABD