ARTERIOCYTE

Helping Patients Heal Faster™

MAGELLAN®

Basic Platelet Separator Kit

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
How Supplied 10

INSTRUCTIONS 10

DISPOSABLES KIT LIMITED WARRANTY – ARTERIOCYTE

MEDICAL SYSTEMS, INC. (OUTSIDE U.S.) 17

MAGELLAN®**Basic Platelet Separator Kit****INTRODUCTION****Indications for Use**

MAGELLAN®¹ Basic Platelet Separator Kit is intended for use only with the MAGELLAN® Autologous Platelet Separator Instrument. See the MAGELLAN® Autologous Platelet Separator Indications for Use.

Contraindications

Use of the MAGELLAN® Basic Platelet Separator Kit is contraindicated for a hemodynamically unstable or hypercoagulable patient.

Use of this product for pediatric patients should be approached carefully. Withdrawing blood from pediatric patients should be at a physician's specific direction, avoiding any significant reduction in the patient's blood volume.

Caution: Medications that adversely affect a patient's coagulation system may inhibit the use of platelet separation system therapy.

Warning: Reprocessing may compromise the structural integrity of the device and/or lead to device failure. Reuse of this single patient use device creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.

Precautions

1. Only Manufacturer-sterilized disposable kits are approved for use with MAGELLAN® Autologous Platelet Separator Instrument. It is important that aseptic technique be used to minimize the possibility of contamination of the disposable components and/or the patient.
2. Store all disposable components in a dry place away from extremes of environmental conditions.
3. Materials used in the MAGELLAN® Basic Platelet Separator Kit may be sensitive to chemicals (such as solvents and

MAGELLAN® is a registered trademark of Arteriocyte Medical Systems, Inc.

certain detergents). Under certain adverse conditions, exposure to these chemicals (including vapors) may cause the plastics to fail or malfunction. Visually inspect the contents of the disposable kit. Should any evidence of damage to the components be found during inspection or setup, do not use the disposable components.

System Description

The MAGELLAN® Basic Platelet Separator Kit consists of sterile fluid pathway single-patient-use components necessary for each platelet separation procedure with the MAGELLAN® Autologous Platelet Separator Instrument.

Cautions

- **! USA Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Refer to the System Operator's Manual supplied with the instrument before performing a platelet separation procedure using the components of this kit.
- Treat all blood and fluids using Universal Precautions for bloodborne pathogens.
- Each platelet separation procedure using the MAGELLAN® Platelet Separator Instrument requires one MAGELLAN® Basic Platelet Separator Kit. The separation chamber and associated tubing can be used with the same patient for up to three complete separation cycles.
- Use only disposable accessories supplied by the Distributor.
- Inspect the kit prior to use. Do not use the kit if any component or the tray is damaged or open.
- Throughout the procedure, make certain all tubing is free of any kinks, twists, or flat areas.
- All components of this kit are sterile fluid pathway single-patient-use. Do not resterilize.

How Supplied

Contents of the MAGELLAN® Basic Platelet Separator Kit (see Figure 1):

- a. One (1) 10-mL BD™² syringe (Syringe 1)
- b. One (1) separation chamber with attached tubing
- c. One (1) 60-mL BD™ syringe (Syringe 2)

Note: If more than one separation cycle will be performed with the same patient, additional 10-mL and 60-mL BD™ syringes are required and can be purchased from the Distributor.


Figure 1. Contents of the MAGELLAN® Basic Platelet Separator Kit.

INSTRUCTIONS

Caution: All components of this kit are sterile fluid pathway single-patient-use. Do not resterilize.

1. Remove the cover from the separator disposables kit tray.
2. Be certain that all components are present and undamaged.

Caution: Do not use the kit if any component or the tray is damaged or open.

3. Remove the separation chamber package and syringes from the carton.
4. Peel open the lid on the chamber package.
5. Holding the platelet separation chamber with the vent facing upward, thread the attached tubing through the center of the chamber caddy (see Figure 2a and b).

² BD™ is a trademark of Becton, Dickinson and Company.


Figure 2. Thread the separation chamber tubing through the center of the chamber caddy.

6. Install the platelet separation chamber into the centrifuge caddy (verify audible snap), make certain that both ends of the chamber are properly located in the caddy notches (see Figure 3a). Place the T-connector into the slot (see Figure 3b), and tubing under the retainer on top of the caddy surface (see Figure 3c). Press the tubing down into the groove on the outer edge of the chamber caddy (see Figure 3d).

Caution: Failure to install the separation chamber properly may result in error codes.


Figure 3. Install the separation chamber into the caddy and place the tubing into the slot and under the retainer.

-
7. Rotate the tubing collar so that its shape aligns with the opening in the support arm. Slide the tubing collar into the support arm (see Figure 4a) and close the latch (see Figure 4b).


Figure 4. Attach chamber tubing collar to the instrument tubing support arm.

8. Press the tubing down into the groove on the support arm and through the notch in the centrifuge ridge (see Figure 5).

Caution: Make certain all tubing is free of any kinks, twists, or flat areas.


Figure 5. Place the tubing into the groove and through the notch in centrifuge ridge.

9. Prepare the syringe pumps on the front of the instrument by rotating the plunger drivers to the open position (see Figure 6a), while sliding them to the lowest points (see Figure 6b).


Figure 6. Within the syringe receptacles, slide the open plunger drivers to the lowest points.

10. Using aseptic technique, draw the appropriate volume of anticoagulant into Syringe 2 (60-mL syringe). Refer to Table 1 for the appropriate volumes of ACD-A and blood.

Table 1. Volumes of ACD-A and corresponding volumes of blood to achieve anticoagulated blood containing ~ 7 parts blood : 1 part ACD-A.

Total Volume of Anticoagulated Blood (mL)	Volume of ACD-A (mL)	Volume of Blood Drawn (mL)
30	4.0	26.0
40	5.0	35.0
50	6.0	44.0
60	8.0	52.0

-
11. Attach Syringe 2 to the blood access site and slowly draw the appropriate volume of blood. Gently mix with ACD-A for thorough distribution. Refer to Table 1 for the appropriate volumes of ACD-A and blood.
Note: Do not exceed 60 mL total volume.
Note: If multiple separation cycles are to be performed with the same patient, repeat steps 10 and 11 to prepare and fill syringes.
 12. Disconnect Syringe 2 from the blood access site.
 13. Remove the luer connector cap from the end of the shorter length of chamber tubing and attach Syringe 1 (10-mL syringe, see Figure 7a).
 14. Evacuate all air from Syringe 1 into the chamber tubing line.
 15. Place Syringe 1 into the appropriate syringe pump receptacle.
 16. Rotate and slide the plunger driver to engage the Syringe 1 plunger (see Figure 7b).
 17. Ensure whole blood is thoroughly mixed using gentle agitation prior to placing in the syringe pump receptacle.
 18. Remove the luer connector cap from the longer length of chamber tubing and attach Syringe 2 (see Figure 7c).
 19. Evacuate all air from Syringe 2 into the chamber tubing line.
Note: Failure to remove air from syringes may compromise the quality of the product.
 20. Place Syringe 2 into the appropriate syringe pump receptacle.

-
21. Rotate and slide the plunger driver to engage the Syringe 2 plunger (see Figure 7d).


Figure 7. Attach tubing, load syringes, and engage syringe plunger drivers.

Caution: Make certain all tubing is free of any kinks, twists, or flat areas.

-
22. Gently stretch the tubing attached to Syringe 2 and slide it completely into the pinch valve opening (see Figure 8).


Figure 8. Stretch Syringe 2 tubing and slide into pinch valve.

Refer to the operator's manual provided with the MAGELLAN® Platelet Separator Instrument for further information on the platelet separation procedure.

DISPOSABLES LIMITED WARRANTY – ARTERIOCYTE MEDICAL SYSTEMS, INC. (OUTSIDE U.S.)¹

THE FOLLOWING LIMITED WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES.

- A. This LIMITED WARRANTY provides assurance for the customer who receives an Arterioocyte Medical Systems® (“AMS”) **MAGELLAN® Basic Platelet Separator Kit** (“Product”), that should the Product fail to function to specification, AMS will issue a credit equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any AMS replacement Product used for that customer. THE WARNINGS CONTAINED IN THE PRODUCT LABELLING ARE CONSIDERED AN INTEGRAL PART OF THIS LIMITED WARRANTY. CONTACT YOUR LOCAL AMS REPRESENTATIVE TO OBTAIN INFORMATION ON HOW TO PROCESS A CLAIM UNDER THIS WARRANTY.
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
- (1) The Product must be used prior to its “Use By” date.
 - (2) The Product must be returned to AMS within sixty (60) days after use and shall be the property of AMS.
 - (3) The Product must not have been altered or subjected to misuse, abuse or accident.
 - (4) The Product must not have been used more than one time for any customer.
 - (5) The Product must be used in conformity with the user instructions in the “User Manual” delivered to the customer together with the Product, of which this LIMITED WARRANTY is an integral part.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
- (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper use or material alteration of the replaced Product.

¹This Limited Warranty is provided by the Manufacturer, Arterioocyte Medical Systems, Inc., 45 South Street, Hopkinton, MA, 01748, USA. It applies only outside of the United States. Areas outside the United States should contact their local Arterioocyte Medical Systems representative for exact terms of the Limited Warranty.

-
- (2) AMS is not responsible for any incidental or consequential damages based on any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort, patent infringement or otherwise.
- D. This LIMITED WARRANTY does not cover those parts that, by their very nature, are likely to deteriorate or which AMS considers should be periodically replaced consistently with normal routine or preventative maintenance requirements.
- E. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part to term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.
- F. No representative, agent, dealer, retailer, or intermediary of AMS shall have authorization to amend the contents of this LIMITED WARRANTY.
- G. The validity, interpretation and performance of the agreement for which this LIMITED WARRANTY is issued, as well as any disputes relating or connected thereto is controlled by and construed under the laws of the State of Ohio, USA.

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