

PRECISE

PRECISION TARGETING FOR PATIENT CARE

Flexibility for different treatment approaches has never been more important. For optimal results in healing, a product solution must be versatile. If procedures don't all fit within the same mold, why should PRP?

Enter Precise™ - the newest advancement in customizable PRP.

Designed to keep patient and physician needs front and center, Precise™ allows for the removal of specific cell types without compromising platelet concentrations.

Whether you're looking for a Leukocyte Rich (LR) or Leukocyte Poor (LP) PRP, Precise™ has you covered.



CONCENTRATE WITHOUT COMPROMISE

VERSATILITY FOR TARGETED CELL FRACTIONS

The system's CollectFlex™ Technology allows for tailored cell compositions with each spin, allowing the system to provide a targeted range of offerings that have never been available before.

Unlike other systems on the market that compromise platelet yields to obtain LP PRP, **Precise can achieve high platelet concentrations** and desired monocyte levels while maintaining **maximum reduction in red blood cells and granulocytes**. Final outputs maintain platelet concentrations above established therapeutic levels².

MULTIPLE OPTIONS MEET THE NEEDS OF TARGET MARKET AND APPLICATION

60mL WB to 7mL PRP	Protocol A - Emcyte Supraphysiologic	Precise Diamond Leukocyte Poor	Precise Sapphire Leukocyte Rich
Platelet Concentration X Baseline	Up to 7.3 xBL	Up to 7.5 xBL	Up to 7.5 xBL
Mononucleated Cells X Baseline	Up to 1.7 xBL	Up to 7 xBL	Up to 8.5 xBL
Granulocytes X Baseline	~ 0.2 xBL	~ 0.2 xBL	~ 0.4 xBL
Hematocrit	~1%	~1%	~2%

STRESS-FREE **SIMPLICITY**

In a single spin, **Precise delivers what you want and removes what you don't** - all while eliminating the interruption of mid-process manipulation.

VALUE IN **FLEXIBILITY**

A Precise disposable can process between 30-60 mL of whole blood in a single, 12-minute cycle. Processing high volumes in a single spin saves time while offering the potential to deliver higher output volumes of PRP.

USER-FRIENDLY SYSTEM

Built with you in mind, the Precise system was designed to decrease the risk of user error. Each single-use, disposable kit is compatible with a thoughtfully designed Precise extraction stand that allows for obtaining finely tuned final volumes of cellular fractions.



PRECISE KIT	PRODUCT CODE
Leukocyte Poor PRP Precise™ Diamond - Cell Concentration Disposable	PRC-DIA-060 PRC-DIA-120
Leukocyte Rich PRP Precise™ Sapphire - Cell Concentration Disposable	PRC-SAP-060 PRC-SAP-120

BRING A TARGETED APPROACH TO CARE
TO YOUR NEXT CASE WITH PRECISE.

TO SPEAK WITH AN ISTO REP,
PLEASE VISIT [ISTOBIO.CO/CONTACT](https://istobio.co/contact)
OR CALL 888.705.ISTO (4786)



References: 1. Data on file, Internal Testing 2. Mautner K et al. A Call for a Standard Classification System for Future Biologic Research: The Rationale for New PRP Nomenclature. PM&R Journal 2015;53-59.
45 South Street • Hopkinton, MA 01748
©2023 Isto Biologics. All rights reserved. Precise™ is a registered trademark of Isto Technologies II LLC. Isto Technologies II LLC and Arterocyte Medical Systems, Inc. are jointly DBA Isto Biologics

istobiologics.com

AMS98306LBL Rev 1