

LIGHTNING™ NEEDLE

**BONE MARROW ASPIRATION
REFINED BY POWER**

DRIVEN BY POWER.

Made for use with any orthopedic drill, Lightning™ Needle assures smooth entry through the strong, bony cortex for quick and easy bone marrow harvesting.

DESIGNED FOR EASE.

The Lightning™ Needle set was designed with performance and ease of use in mind. While new designs of the hub and trocar tip aid in draw efficiency, large bevel-edged fenestrations allow for bone marrow to gently flow through the cannula in an effort to increase flow and decrease pressure on cells.



Cannula + Trocar



Plunger



Adapter



Eight bevel-edged fenestrations allow greater space for marrow flow.



Extended trocar tip ensures that the cutting edge is always engaged.

Grooves on the cannula and trocar hubs increase grip potential.